The human epithelial lung cell line "A549" as a suitable tool in *in vitro* inhalation toxicity testing



Detlef Ritter, Jan Knebel Fraunhofer ITEM, Germany

Background

- Air-lifted interface (ALI) culture methods represent a state-of-the-art approach for investigations on biological effects of inhalable, gases, vapors or aerosols.
- Usually, such investigations aim at a high in vivo relevance by application of more complex cell models such as primary cells, 3D co-cultures or similar.
- Due to a demand for high numbers of cultures for routine testing especially in toxicity screening studies, "cheap and easy-to-use" cell culture models are urgently needed.

Objectives

- The human alveolar epithelial type-II like cell line A549 represents a comprehensively characterized culture model (> 20.000 references)
- Although lacking tight-junctions and also with a carcinogenic origin, many studies have also confirmed primary characteristics (CYP-expression, surfactant biosynthesis, etc.)
- The objective of this overview is to highlight the qualities of A549 ALI cultures in inhalation toxicity testing *in vitro* with respect to applicability, in vivo relevance and reproducibility
- Summary of selected inhalation in vitro data from studies of the ITEMs lab during recent years

Inhalation toxicity testing in vitro: Methods



Applicability to gases, vapors and aerosols from test items, experimental or environmental atmospheres



Dry Particle aerosols *Testing of fungicide aerosols*



Experimental or environmental aerosols and atmospheres *Testing of combustion products: battery failure or tobacco smoke*





In vivo relevance by *in vitro* ↔ *in vivo* correlations

Approach a) ...by referencing to relevant positive controls Human inhalation irritants



120

Dry particle aerosol inhalation Gas/vapor inhalation *in vivo / in vitro data* in vivo / in vitro data rats) y = 0,2249x^{0,4615} $R^2 = 0,9207$ **in vivo** inhalation, ¹ C50 [mg/l] $y = 2,1988\ln(x) - 4,8369$ $R^2 = 0,9968$ vivo i LC Chlorothalonil Ŀ, Mancozeb A (commercial fungicide) 0 B (commercial fungicide) •••••• Fit (all data) 100 10000 100000 1000 1000000 in vitro ⁵⁰ Cell based in vitro Toxicity ED₅₀ [ppm x h x ml/min x cm⁻¹] ("in vitro inhalation", human lung cells) EC50 [μg/cm²] ▲ tert buty hydroperoxide ● isobutanol

Reproducibility of *in vitro* inhalation toxicity data over years

Gas/vapor exposure dose-response data

Dry particle aerosol exposure EC_{50} data

20

styrene

▲ toluene

100



Approach b) ...by referencing to in vivo inhalation data





Conclusions

- A549 human lung ALI cultures represent a relatively easy-to-use biological model that can be made available in the lab in large culture numbers (e.g. > 100/week) at moderate costs.
- Inhalation toxicity studies using A549 human lung ALI cultures have been shown to lead to a relevant prediction of in vivo toxicity.
- If applied under appropriate conditions, considerable reproducibility can be achieved, necessary for establishment of *in vitro* toxicity databanks and historical lab controls.
- These characteristics qualify A549 ALI cultures as biological models in first stages of tiered in vitro inhalation approaches, mechanistic studies may follow applying more complex cell models.

detlef.ritter@item.fraunhofer.de www.item.fraunhofer.de www.prit-systems.de